

## Product Datasheet

### CD11a antibody (FITC) (orb623056)

|                            |                                                                                             |
|----------------------------|---------------------------------------------------------------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody against CD11a conjugated to FITC                                  |
| <b>Species/Host</b>        | Mouse                                                                                       |
| <b>Reactivity</b>          | Human                                                                                       |
| <b>Conjugation</b>         | FITC                                                                                        |
| <b>Tested Applications</b> | FC                                                                                          |
| <b>Immunogen</b>           | Leucocytes from a patient suffering from a LGL-type leukemia                                |
| <b>Target</b>              | CD11a                                                                                       |
| <b>Preservatives</b>       | 0.09% Sodium Azide (NaN <sub>3</sub> ) 0.2% Bovine Serum Albumin. Phosphate buffered saline |
| <b>Form/Appearance</b>     | Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid              |
| <b>Storage</b>             | +4°C DO NOT FREEZE                                                                          |
| <b>Note</b>                | For research use only                                                                       |
| <b>Isotype</b>             | IgG1                                                                                        |
| <b>Clonality</b>           | Monoclonal                                                                                  |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                             |