

Product Datasheet

CD11a antibody (RPE) (orb623057)

Description	Mouse monoclonal antibody against CD11a conjugated to RPE
Species/Host	Mouse
Reactivity	Human
Conjugation	RPE
Tested Applications	FC
Immunogen	Leucocytes from a patient suffering from a LGL-type leukemia
Target	CD11a
Preservatives	0.1% Sodium Azide (NaN ₃) 0.2% Bovine Serum Albumin. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to R. Phycoerythrin (RPE) - liquid
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.