

Product Datasheet

CD195 antibody (RPE) (orb623089)

Description	Mouse monoclonal antibody against CD195 conjugated to RPE
Species/Host	Mouse
Reactivity	Human
Conjugation	RPE
Tested Applications	FC
Immunogen	CCR5 peptide (Met1-Lys22) KLH conjugate
Target	CD195
Preservatives	0.1% Sodium Azide (NaN ₃) 0.2% Bovine Serum Albumin. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to R. Phycoerythrin (RPE) - liquid
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.