

Product Datasheet

Human G-CSF protein (orb623133)

Description	Recombinant human G-CSF protein
Reactivity	Human
Tested Applications	FA
Target	G-CSF
Preservatives	1.0% Trehalose. 20 mM Phosphate Buffer, 0.1 M Sodium Chloride
Form/Appearance	Purified recombinant protein - lyophilized
Concentration	500 µg/ml after reconstitution
Storage	-20°C only (ship frozen)
Note	For research use only
Clonality	Recombinant
Purity	Rec. Protein
Expiration Date	6 months from date of receipt.