

## Product Datasheet

### Mouse EGF protein (orb623136)

|                            |                                            |
|----------------------------|--------------------------------------------|
| <b>Description</b>         | Recombinant mouse EGF protein              |
| <b>Reactivity</b>          | Mouse                                      |
| <b>Tested Applications</b> | FA, WB                                     |
| <b>Target</b>              | EGF                                        |
| <b>Preservatives</b>       | None present                               |
| <b>Form/Appearance</b>     | Purified recombinant protein – lyophilized |
| <b>Concentration</b>       | 1.0 mg/ml after reconstitution             |
| <b>Storage</b>             | -20°C only (ship +4°C)                     |
| <b>Note</b>                | For research use only                      |
| <b>Clonality</b>           | Recombinant                                |
| <b>Purity</b>              | Rec. Protein                               |
| <b>Expiration Date</b>     | 6 months from date of receipt.             |