

Product Datasheet

Goat anti Rat IgG/IgM (H/L) antibody (FITC) (orb623140)

Description	Goat polyclonal anti Rat IgG/IgM (H/L) conjugated to FITC
Species/Host	Goat
Reactivity	Rat
Conjugation	FITC
Tested Applications	FC
Immunogen	Purified rat IgM and IgG
Target	IgG IgM
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Concentration	IgG concentration 1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.