

Product Datasheet

Goat anti Rabbit IgG/IgM antibody (RPE) (orb623149)

Description	Goat polyclonal anti Rabbit IgG/IgM conjugated to RPE
Species/Host	Goat
Reactivity	Rabbit
Conjugation	RPE
Tested Applications	FC
Immunogen	Normal rabbit IgG + IgM.
Target	IgG IgM
Preservatives	0.1% Sodium Azide (NaN ₃)Stabilizing agent (sucrose). Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to R. Phycoerythrin (RPE) - liquid
Concentration	IgG concentration 0.5 mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.