

Product Datasheet

Parvalbumin Antibody (orb623238)

Description	Parvalbumin Antibody
Species/Host	Gallus
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Full-length recombinant human Parvalbumin Protein
Form/Appearance	Affinity Purified
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	Chicken IgY
Uniprot ID	P20472
Entrez	5816
Expiration Date	12 months from date of receipt.