

Product Datasheet

Nestin Antibody (orb623241)

Description	Nestin Antibody
Species/Host	Gallus
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Recombinant construct, amino acids 317-630 of the human protein expressed in E coli.
Storage	Store at 4°C short term. Store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	Chicken IgY
Uniprot ID	P48681
Entrez	10763
Dilution Range	IF/IHC: 1:500-1:1.000 Western Blot: 1:5,000
Expiration Date	12 months from date of receipt.