

Product Datasheet

mCherry - Mouse Antibody (orb623264)

Description	mCherry - Mouse Antibody
Species/Host	Mouse
Reactivity	Mammal
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Antibody was raised in mouse against recombinant full length His tagged mCherry purified from E. coli.
Form/Appearance	Affinity Purified
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	Mouse IgG
Expiration Date	12 months from date of receipt.