

Product Datasheet

Aurora B Kinase Antibody (orb623268)

Description	Aurora B Kinase Antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Full length recombinant human aurora B expressed in and purified from E. coli.
Form/Appearance	Affinity Purified
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	Mouse IgG
Uniprot ID	Q96GD4
Entrez	9212
Expiration Date	12 months from date of receipt.