

Product Datasheet

Green Fluorescent Protein Antibody (orb623283)

Description	Green Fluorescent Protein Antibody
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant Aequoria coerulescens GFP purified from E. coli
Form/Appearance	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	Mouse IgM
Uniprot ID	P42212
Entrez	7011691
Dilution Range	IF/ICC: 1:1,000-5,000 Immunohistochemistry: 1:1,000-5,000 Western Blot: 1:1,000-5,000
Expiration Date	12 months from date of receipt.