

## Product Datasheet

# Synaptosomal Associated Protein 25/SNAP25 peptide (orb623293)

|                        |                                                                                          |
|------------------------|------------------------------------------------------------------------------------------|
| <b>Description</b>     | Synaptosomal Associated Protein 25/SNAP25 peptide                                        |
| <b>Form/Appearance</b> | Frozen                                                                                   |
| <b>Storage</b>         | Store at 4°C short term. Store at -20°C for maximum stability. Avoid freeze-thaw cycles. |
| <b>Note</b>            | For research use only                                                                    |
| <b>Expiration Date</b> | 6 months from date of receipt.                                                           |