

Product Datasheet

Human Bromodomain Containing 1 protein (orb623405)

Description	Human Bromodomain Containing 1 protein
Reactivity	Human
Form/Appearance	Liquid
Storage	Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.
Note	For research use only
Uniprot ID	095696
Entrez	23774
Expiration Date	12 months from date of receipt.