

Product Datasheet

Human Breast Cancer Metastasis Suppressor 1 protein (orb623408)

Description	Human Breast Cancer Metastasis Suppressor 1 protein
Reactivity	Human
Form/Appearance	Liquid
Storage	Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.
Note	For research use only
Uniprot ID	Q9HCU9
Entrez	25855
Expiration Date	12 months from date of receipt.