

Product Datasheet

Neurofilament NF-M Antibody (orb623504)

Description	Neurofilament NF-M Antibody
Species/Host	Rabbit
Reactivity	Feline, Gallus, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IHC, WB
Immunogen	Recombinant fusion protein containing the extreme C-terminal segment of rat NF-M
Form/Appearance	Whole serum
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Uniprot ID	P07197
Entrez	4741
Dilution Range	ICC 1:100; IHC 1:100; WB 1:1,000
Expiration Date	12 months from date of receipt.