

Product Datasheet

Laminin 111 Antibody (orb623506)

Description	Laminin 111 Antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Uniprot ID	P25391
Entrez	284217
Expiration Date	12 months from date of receipt.