

Product Datasheet

Nestin Antibody (orb623510)

Description	Nestin Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Amino acids 315-630 of the human protein
Form/Appearance	Whole serum
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Uniprot ID	P48681
Entrez	10763
Expiration Date	12 months from date of receipt.