

Product Datasheet

IBA1 antibody (orb623513)

Description	IBA1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Peptide identical to part of the C-terminal of human IBA1 coupled to KLH.
Form/Appearance	Whole serum
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
Uniprot ID	P55008
Entrez	199
Dilution Range	Western Blot: 1:1,000-5,000. Immunocytochemistry/Immunofluorescence: 1:2,000-5,000
Expiration Date	12 months from date of receipt.