

Product Datasheet

Ankyrin 3 Antibody (orb623516)

Description	Ankyrin 3 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	C-terminal 398 amino acids of human ankyrin 3 expressed in and purified in E. coli
Form/Appearance	Liquid
Storage	Store at 4°C short term. Store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Uniprot ID	Q12955
Entrez	288
Dilution Range	IF/ICC: 1:2,000-1:5,000 IHC: 1:2,000-1:5,000 WB: 1:1,000-5,000
Expiration Date	12 months from date of receipt.