

Product Datasheet

Calbindin D-28K Antibody (orb623520)

Description	Calbindin D-28K Antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC
Immunogen	Calbindin D -28K purified from bovine cerebellum
Form/Appearance	Whole serum
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Uniprot ID	P05937
Entrez	793
Expiration Date	12 months from date of receipt.