

Product Datasheet

Parvalbumin Antibody (orb623521)

Description	Parvalbumin Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Primate, Rat
Conjugation	Unconjugated
Tested Applications	IHC
Form/Appearance	Whole serum
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Uniprot ID	P20472
Entrez	5816
Expiration Date	12 months from date of receipt.