

Product Datasheet

CD66e antibody (APC) (orb623577)

Description	Mouse monoclonal antibody against CD66e conjugated to APC
Reactivity	Human
Conjugation	APC
Tested Applications	FC
Immunogen	Human carcinoembryonic antigen (CEA; CEACAM5)
Target	CD66e
Preservatives	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	0.1 mg/ml
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified antibody is conjugated with activated allophycocyanin (APC) under optimum conditions and unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.
Uniprot ID	P06731
Entrez	1048
Expiration Date	12 months from date of receipt.