


Product Datasheet


Syndecan 1 antibody (APC) (orb623578)

Description Mouse monoclonal antibody against Syndecan 1


Reactivity	Human, Primate, Rat
Conjugation	APC
Tested Applications	FC
Immunogen	A mixture of U266 and XG-1 human myeloma cell lines
Target	CD138
Preservatives	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: The reagent is designed for analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 106 cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests.
Isotype	Mouse IgG1 kappa
Clonality	Monoclonal
Purity	Purified antibody is conjugated with activated allophycocyanin (APC) under optimum conditions and unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.
Uniprot ID	P18827
Entrez	6382
Dilution Range	Flow cytometry: The reagent is designed for analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 106 cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests.
Expiration Date	12 months from date of receipt.



Staining pattern of Anti-human CD38 FITC...



Staining pattern of Anti-human CD38 FITC...



Reactivity of Anti-human CD38 FITC antib...