



Product Datasheet

IgE antibody (Biotin) (orb623580)

Description	Mouse monoclonal antibody against IgE conjugated to Biotin
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA, FC
Immunogen	Purified human IgE.
Target	IgE
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-12 µg/ml
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified antibody is conjugated with biotin LC-NHS ester under optimum conditions and unconjugated antibody and free biotin are removed by size-exclusion chromatography.
Dilution Range	Flow cytometry: Recommended dilution: 1-12 µg/ml
Expiration Date	12 months from date of receipt.



Separation of human IgE positive basophi...



Flow cytometry surface staining pattern ...