



Product Datasheet

CD16 antibody (PE) (orb623585)

Description	Rat monoclonal antibody against CD16 conjugated to PE
Reactivity	Mouse
Conjugation	PE
Tested Applications	FC
Immunogen	Murine pre-B cells
Target	CD16/CD32
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	0.5 mg/ml
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 2 µg/ml.
Isotype	Rat IgG2a lambda
Clonality	Monoclonal
Purity	Purified antibody is conjugated with R-phycoerythrin (PE) under optimum conditions. Unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.
Dilution Range	FACS: 2 µg/ml
Expiration Date	12 months from date of receipt.



Separation of murine CD16/32 positive ce...



Flow cytometry surface staining pattern ...