


Product Datasheet

CD21 antibody (PerCP) (orb623596)


Description

Mouse monoclonal antibody against CD21 conjugated to

Reactivity	Bovine, Canine, Human, Porcine
Conjugation	PerCP
Tested Applications	FC
Immunogen	IM9 human B-lymphoblastoid cell line
Target	CD21
Preservatives	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: The reagent is designed for analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 106 cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests.
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified antibody is conjugated with activated Peridinin-Chlorophyll Protein (PerCP) under optimum conditions and unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.
Uniprot ID	P20023
Entrez	1380
Dilution Range	Flow cytometry: The reagent is designed for analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 106 cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests.
Expiration Date	12 months from date of receipt.



Separation of human CD21 lymphocytes (re...



Flow cytometry surface staining pattern ...