


Product Datasheet

Rabbit IGF-1 ELISA Kit (orb623709)

Description	Rabbit IGF-1 ELISA Kit
Reactivity	Rabbit
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: G49-A118
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
MW	21841 MW
Uniprot ID	P05019
Sensitivity	50 pg/ml
Expiration Date	Please enquire.



Rabbit IGF-1
PicoKine
ELISA Kit
standard...