

Product Datasheet

SAE2/UBA2 Antibody (monoclonal, 5B13) (orb623715)

Description

SAE2/UBA2 Antibody (monoclonal, 5B13)

Species/Host

Mouse

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

IHC, WB

Applications

Immunogen

E. coli-derived human SAE2/UBA2 recombinant protein (Position: E449-K564).

Form/Appearance

Lyophilized

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Add 0.2ml of distilled water will yield a concentration of 500µg/ml.

Isotype

Mouse IgG2b

Clonality

Monoclonal

MW

90 kDa

Uniprot ID


[Q9UBT2](#)

Dilution Range


Western blot, 0.1-0.5µg/ml, Human, Mouse, Rat
Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml, Human, By Heat

Expiration Date


12 months from date of receipt.



IHC analysis of SAE2/UBA2 using anti-SAE...



IHC analysis of SAE2/UBA2 using anti-SAE...



IHC analysis of SAE2/UBA2 using anti-SAE...