



Product Datasheet

IRS1 (phospho-Ser789) antibody (orb6238)


Description	IRS1 (phospho-Ser789) antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC-P
Immunogen	KLH conjugated synthesised phosphopeptide derived from human IRS1 around the phosphorylation site of Ser789 (SS(p-S)SG)
Target	phospho-IRS1 (Ser789)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	132 kDa
Uniprot ID	P35568
Dilution Range	IHC-P=1:100-500, ICC=1:100, Flow-Cyt=1µg /test, IHC-F=1:100-500, IF=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



Flow cytometry of HeLa(blue line).Primar...



ICC analysis of A549 cell using phospho-...



IHC-P analysis of Rat stomach using p-IR...