


Product Datasheet

Synaptophysin/SYP Antibody (orb623902)


Description

Synaptophysin/SYP Antibody


Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	A synthetic peptide corresponding to a sequence of human Synaptophysin/SYP (DLSIEVEFEYPFRLHQVYFDAPTCR).
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Add 0.2ml of distilled water will yield a concentration of 500ug/ml.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	38 kDa
Uniprot ID	P08247
Dilution Range	Western blot, 0.1-0.25µg/ml, Mouse, Rat Immunohistochemistry(Paraffin-embedded Section, 2-5 µg/ml, Human, Mouse, Rat Immunocytochemistry/Immunofluorescence, 5µg/ml, Human Immunofluorescence, 5µg/ml, Mouse, Rat
Expiration Date	12 months from date of receipt.



Western blot analysis of SYP using anti-...



IF analysis of SYP using anti-SYP antibody...



IHC analysis of SYP using anti-SYP antibody...