

Product Datasheet

IRS1 (phospho-Tyr895) antibody (orb6240)

Description	IRS1 (phospho-Tyr895) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthesised phosphopeptide derived from human IRS1 around the phosphorylation site of Tyr895 (GE(p-Y)VN)
Target	phospho-IRS1 (Tyr895)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	132 kDa
Uniprot ID	P35568
Expiration Date	12 months from date of receipt.