

Product Datasheet

Vimentin antibody (orb624175)

Description	Mouse monoclonal antibody to Vimentin
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-Fr, IHC-P, WB
Immunogen	Please contact us for more information.
Preservatives	Culture medium with 0.01M Sodium Azide
Concentration	30 µg/ml
Storage	2-8°C
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	IHC-P (1:100),IHC-fr (1:100),,
Expiration Date	12 months from date of receipt.