

Product Datasheet

Keratin 4 antibody (orb624177)

Description	Mouse monoclonal antibody to Keratin 4
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-Fr, WB
Immunogen	Please contact us for more information.
Preservatives	Culture medium with 0.01M Sodium Azide
Concentration	15 µg/ml
Storage	2-8°C
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	IHC-Fr (1:20),,
Expiration Date	12 months from date of receipt.