

Product Datasheet

Mouse Epidermal Growth Factor antibody (orb624193)

Description	Mouse monoclonal antibody to Mouse Epidermal Growth Factor
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-Fr, IHC-P, WB
Immunogen	Please contact us for more information.
Preservatives	Culture medium with 0.01M Sodium Azide
Concentration	20 µg/ml
Storage	2-8°C
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	IHC-P (1:5), IHC-fr (1:5),,
Expiration Date	12 months from date of receipt.