

Product Datasheet

MRP1 antibody (orb624201)

Description	Mouse monoclonal antibody to MRP1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC-Fr, IHC-P, WB
Immunogen	Please contact us for more information.
Preservatives	Serum free tissue culture supernatant with 0.7% rAlbumin and 0.1% sodium azide
Concentration	250 µg/ml
Storage	2-8°C
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	,,,IHC-P (1:20-1:50)
Expiration Date	12 months from date of receipt.