

Product Datasheet

**Monkey Serum Albumin protein [Out of stock]
(orb624208)**

Description

serum albumin (SA) is also known as ALB, which is

kDa	M	R
116.0		

Reactivity

Monkey

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Tested

Applications

SDS-PAGE

Target

Serum Albumin

Preservatives

PBS, pH 7.4

Form/Appearance

Powder

Storage

-20°C

Tag

C-10×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 67.9 kDa. The protein migrates as 66 kDa under reducing (R) condition (SDS-PAGE).

Protein Sequence

NP_001306284.1

Purity

95%

Source

Cynomolgus Serum Albumin, His Tag (CSA-C52H4) is expressed from human 293 cells (HEK293). It contains AA Asp 25 - Ala 608 (Accession # A2V9Z4-1).

MW


67.9 kDa

NCBI

[NP_001306284.1](https://www.ncbi.nlm.nih.gov/nuccore/NP_001306284.1)

Expiration Date

6 months from date of receipt.



SDS-Page
analysis of
Monkey
Serum
Albumi...