

Product Datasheet

**Rabbit Serum Albumin protein [Out of stock]
(orb624209)**

Description

serum albumin (SA) is also known as

plasma and has a good binding ...

Reactivity

Rabbit

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Tested

SDS-PAGE

Applications
Target

Serum Albumin

Preservatives

PBS, pH 7.4

Form/Appearance

Powder

Storage

-20°C

Tag

C-10×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 67.9 kDa. The protein migrates as 60-70 kDa under reducing (R) condition (SDS-PAGE).

Protein Sequence

NP_001075813.1

Purity

95%

Source

Rabbit Serum Albumin, His Tag (RSA-R52H6) is expressed from human 293 cells (HEK293). It contains AA Glu 25 - Gly 608 (Accession # P49065-1).

MW

67.9 kDa

NCBI
[NP_001075813.1](#)
Expiration Date

6 months from date of receipt.