


Product Datasheet

Human PGP9.5 ELISA kit (orb624277)

Description	Human PGP9.5(Ubiquitin carboxyl-terminal hydrolase isozyme L1) ELISA Kit
Reactivity	Human
Range	78-5000pg/ml
Target	UCHL1
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P09936
Sensitivity	46.9pg/ml
Expiration Date	Please enquire.



Human
UCHL1(Ubiquitin
carboxyl-
terminal ...