


## Product Datasheet

### Human sFMC ELISA kit (orb624278)

|                        |                       |
|------------------------|-----------------------|
| <b>Description</b>     | Human sFMC ELISA kit  |
| <b>Reactivity</b>      | Human                 |
| <b>Range</b>           | 1.563-100ug/ml        |
| <b>Target</b>          | sFMC                  |
| <b>Storage</b>         | 2-8 °C for 6 months   |
| <b>Note</b>            | For research use only |
| <b>Sensitivity</b>     | 0.938ug/ml            |
| <b>Expiration Date</b> | Please enquire.       |



Human  
sFMC(Soluble  
Fibrin  
Monomer  
Comple...