


## Product Datasheet

### Human beta 1-ARA ELISA kit (orb624280)

|                        |                            |
|------------------------|----------------------------|
| <b>Description</b>     | Human beta 1-ARA ELISA kit |
| <b>Reactivity</b>      | Human                      |
| <b>Range</b>           | 0.781-50ng/ml              |
| <b>Storage</b>         | 2-8 °C for 6 months        |
| <b>Note</b>            | For research use only      |
| <b>Uniprot ID</b>      | <b>P08588</b>              |
| <b>Sensitivity</b>     | 0.469ng/ml                 |
| <b>Expiration Date</b> | Please enquire.            |



Human  
beta 1-  
ARA(beta  
1  
adrenergic  
recep...