


## Product Datasheet

### 4-HNE ELISA kit (orb624366)

|                        |                       |
|------------------------|-----------------------|
| <b>Description</b>     | 4-HNE ELISA kit       |
| <b>Reactivity</b>      | All                   |
| <b>Range</b>           | 31.25-2000pg/ml       |
| <b>Storage</b>         | 2-8 °C for 6 months   |
| <b>Note</b>            | For research use only |
| <b>Sensitivity</b>     | 18.75pg/ml            |
| <b>Expiration Date</b> | Please enquire.       |



4-HNE (4-Hydroxynonenal)  
ELISA Kit