

Product Datasheet

14-3-3 antibody (orb624373)

Description	Rabbit polyclonal antibody to 14-3-3
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, IP, WB
Immunogen	tyrosine 3-monooxygenase/tryptophan 5-monooxygenase activation protein, theta polypeptide
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	28 kDa
Uniprot ID	P27348
Dilution Range	WB: 1:1000-1:10000, IP: 1:200-1:2000, IHC-P: 1:20-1:200, IF: 1:20-1:200
Expiration Date	12 months from date of receipt.