

Product Datasheet

ABCB9 antibody (orb624398)

Description	Rabbit polyclonal antibody to ABCB9
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	ATP-binding cassette, sub-family B(MDR/TAP), member 9
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	84 kDa
Uniprot ID	Q9NP78
Dilution Range	WB: 1:200-1:2000, IF: 1:10-1:100
Expiration Date	12 months from date of receipt.