

Product Datasheet

ACE antibody (orb624440)

Description	Rabbit polyclonal antibody to ACE
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	angiotensin I converting enzyme(peptidyl-dipeptidase A) 1
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	150 kDa, observed 170 kDa
Uniprot ID	P12821
Dilution Range	WB: 1:500-1:5000, IP: 1:500-1:5000
Expiration Date	12 months from date of receipt.