

Product Datasheet

ACRV1 antibody (orb624460)

Description	Rabbit polyclonal antibody to ACRV1
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IP, WB
Immunogen	acrosomal vesicle protein 1
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	28 kDa
Uniprot ID	P26436
Dilution Range	WB: 1:1000-1:10000, IP: 1:1000-1:10000, IF: 1:20-1:2000
Expiration Date	12 months from date of receipt.