

Product Datasheet

ACSM5 antibody (orb624468)

Description	Rabbit polyclonal antibody to ACSM5
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	acyl-CoA synthetase medium-chain family member 5
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	23 kDa, observed 65 kDa
Uniprot ID	Q6NUN0
Dilution Range	WB: 1:200-1:2000, IHC-P: 1:20-1:200
Expiration Date	12 months from date of receipt.