

Product Datasheet

ACSS1 antibody (orb624469)

Description	Rabbit polyclonal antibody to ACSS1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	acyl-CoA synthetase short-chain family member 1
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	75 kDa
Uniprot ID	Q9NUB1
Dilution Range	WB: 1:1000-1:10000, IP: 1:200-1:2000
Expiration Date	12 months from date of receipt.