

Product Datasheet

ACTR1B antibody (orb624480)

Description	Rabbit polyclonal antibody to ACTR1B
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	ARP1 actin-related protein 1 homolog B, centractin beta(yeast)
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	42 kDa
Uniprot ID	P42025
Dilution Range	WB: 1:200-1:2000, IHC-P: 1:20-1:200, IF: 1:20-1:200
Expiration Date	12 months from date of receipt.