

Product Datasheet

ADARB1 antibody (orb624504)

Description	Rabbit polyclonal antibody to ADARB1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	adenosine deaminase, RNA-specific, B1(RED1 homolog rat)
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	81 kDa
Uniprot ID	P78563
Dilution Range	WB: 1:500-1:5000, IP: 1:500-1:5000
Expiration Date	12 months from date of receipt.